

# VIRAL INACTIVATION FLUID



**tes**corn



## INGREDIENTS

Tris-HCL , EDTA, Protein Denaturant , Lipid Solvent

## PHYSICAL PROPERTIES

Liquid, colorless

## DESCRIPTION

The product is used for the transfer of enveloped and unenveloped viruses in a inactive condition for detection and research purposes.

## PRODUCT SPECIFICATION

The product is used in molecular virology, molecular biology, molecular genetic and chemistry fields. It is suitable for RNA/DNA isolation procedures and it does not contain DNase, RNase, exonucleases and endonucleases.

Only for laboratory use.

## PROCEDURE

1. Immerse nasal/ nasopharyngeal swab sample directly in the inactivation solution and rotate for 3-5 times.
2. Discard the swab as medical waste.
3. Incubate inactivation solution at room temperature for 10 minutes.
4. Virally inactivated sample can be used for DNA/RNA isolation and all following molecular studies including Real Time PCR.

## CAUTIONS

- In case of contact, it may irritate skin.
- Do not use without gloves.
- In case of contact, immediately wash skin with copious amounts of water.
- Do not use expired products.
- Shelf-life of the product is 18 months.
- Please contact us for your questions.

## STORAGE

Room temperature

## DATE OF ISSUE

20.03.2020



### TESSCORN SEASIA PTE. LIMITED

596b Ang Mo Kio Street 52  
#17-327, Singapore (562596)

enquiries@tesscorn-seasia.com

### TESSCORN USA LLC.

3895 Lawson Dr. ,  
Troy, MI 48084

Cell: +1 734 658 0450  
raghu@tesscorn-usa.com

### TESSCORN SYSTEMS INDIA PVT LTD

1285, 5th Main, 17th Cross  
HSR Sector 7, Bangalore 560068, India

Cell: +91 8861751155, +91 9900240316  
info@tesscorn.com